

## SPECTROPHOTOMETER UV-VIS DOUBLE BEAM. Model UV-6900



Spectrophotometer belonging to a new generation of dual-beam instruments with 2 detectors that simultaneously measure the reference and the sample to optimize accuracy, sensitivity, low stray light and a bandwidth of up to 0.5 nm.

A suitable instrument for clinical, pharmaceutical and biochemical laboratories as well as current applications and quantitative analysis, kinetics, scanings at various wavelengths and analysis of DNA / Protein. All operable directly on the instrument and displayed on the LCD screen or transferred to PC, data, and curve.

- ◆ Auto - calibration baseline, wavelength and dark current
- ◆ Equipped with parallel printer port and one USB for PC. Included software
- ◆ Available a variety of accessories to complete the versatility of the spectrophotometer.
- ◆ New layout monochromator diffraction grating with 1200 lines which ensures high durability and stability of the system
- ◆ Software Update by internet
- ◆ Real-time results, date and time
- ◆ 320x240mm LCD big screen

Parameters	Features	Parameters	Features
<b>Optical principle</b>	Double beam	<b>Baseline flatness</b>	±0.001Abs
<b>Wavelength range</b>	190-1100nm	<b>Noise</b>	±0.001 Abs
<b>Wavelength accuracy</b>	±0.3nm	<b>Band width</b>	0.5 / 1.0 / 2.0 / 4, to select
<b>Wavelength repeatability</b>	±0.2nm	<b>Stray light</b>	<0.05%T (220nm, 360nm)
<b>Photometric accuracy</b>	±0.3%T	<b>Out</b>	USB
<b>Measuring ranges</b>	Abs: -0.3 to 3.0 %T: 0 to 200%T Conc: 0 to 9999C	<b>Baseline stability</b>	0.002 Abs/h ( at 500nm, after 2h warm up
<b>Photometric repeatability</b>	±0.2%T	<b>Light source</b>	W1, D2
<b>Instrument control</b>	Autonomous and PC	<b>Detectors</b>	Silicon photodiodes
<b>Weight</b>	30Kg	<b>Printer port</b>	Parallel
<b>Power</b>	220VAC 50/60Hz	<b>Size</b>	625w x 430d x 206h mm

### Accessories



8-position 10mm. automatic cell - holder. Code 6.6900.01



Micro-printer. Code 6.6900.07

## SOFTWARE UV/VIS ANALYST (included). Features

### Multi - Wavelength

Up 32 wavelengths can be selected and multiple samples can be measured

### Quantitative

- 1- Coefficient method
- 2- Standard curve  
Up to 10 standards samples may be used to establish a curve. Four methods for fitting a curve through the calibration points: Linear fit, Linear fit through zero, Square fit and Cubic fit

### Wavelength Scan

- 1- The wavelength scan intervals are: 0.1 - 0.2 - 0.5 - 1 - 2.5
- 2- Scan speed: High - Medium and Low, from 100 to 3000nm/min.

### DNA /Protein Test

Concentration and DNA purity are quickly and easily calculated:

Absorbance ratios 260nm / 280nm with optional subtracted absorbance at 320nm

DNA concentration =  $62.9 \times A_{260} - 36.0 \times A_{280}$

Protein concentration =  $1552 \times A_{260} - 757.3 \times A_{280}$

### Kinetics

Calculate reaction rates. Graphical time / real time absorbance

Total time up to 12 hours with reading selectable interval 0.5 - 1 - 2 - 5 - 10 - 30 seconds and 1 minute

Further manipulation with re-scaling, curve tracking and selection of the part of the curve required for the rate calculation

The speed is calculated using a linear regression algorithm before multiplying by the factor entered

Code	Description
<b>6.6900.00*</b>	Spectrophotometer Double Beam . Model UV 6900 with software, cable USB, 2-cell cuvette 10mm holder, 4 glass cuvettes 10mm, 2 quartz cuvettes 10mm. Band width 0.5 / 1.0 / 2.0 / 4.0 to select
<b>6.6900.01</b>	8-position 10mm. automatic cell - holder
<b>6.6900.03</b>	Halogen lamp
<b>6.6900.08</b>	Deuterium lamp
<b>6.6900.06</b>	Cuvette with lid, quartz, 3.5ml, path light 10mm, 2pcs
<b>6.6900.05</b>	Cuvette with lid, glass, 3.5ml, path light 10mm, 4pcs
<b>8.9749.00*</b>	PC, with monitor TFT 19" for UV6900
<b>6.6900.07*</b>	Micro-printer for dates

Additional accessories on request

Distributor: **Tajhiz Pouyesh Amitis**    [www.tpaco.ir](http://www.tpaco.ir)    6641369