



Instruction Manual

**EC214
EC215 - EC215R
Bench
Conductivity
Meters**



HANNA
instruments
<http://www.hannainst.com>

CE
These Instruments are in
Compliance with the CE Directives

Dear Customer,

Thank you for choosing a Hanna Instruments Product.

Please read this instruction manual carefully before using the instrument.

This manual will provide you with all the necessary information for the correct use of the instrument, as well as a precise idea of its versatility in a wide range of applications.

These instruments are in compliance with **CE** directives.

TABLE OF CONTENTS

PRELIMINARY EXAMINATION	3
GENERAL DESCRIPTION	3
FUNCTIONAL DESCRIPTION EC 214	4
SPECIFICATIONS EC 214	5
FUNCTIONAL DESCRIPTION EC 215	6
SPECIFICATIONS EC 215	7
OPERATIONAL GUIDE	8
CALIBRATION	10
CONDUCTIVITY VS. TEMPERATURE CHART	13
DETERMINATION OF THE TEMPERATURE	
COEFFICIENT OF A SOLUTION (EC 215 & EC 215R)	14
PROBE MAINTENANCE	15
ACCESSORIES	16
WARRANTY	18
CE DECLARATION OF CONFORMITY	19



PRELIMINARY EXAMINATION

Remove the instrument from the packing material and examine it carefully to make sure that no damage has occurred during shipping. If there is any damage, notify your dealer.

Each meter is supplied complete with:

- Conductivity Probe (HI76300 for EC 214, HI76303 for EC 215 and EC 215R);
- Instruction Manual;
- 12VDC adapter (HI710005 or HI710006);
- Dust cover.

Note: Save all packing material until you are sure that the instrument functions correctly. All defective items must be returned in the original packing with the supplied accessories.

GENERAL DESCRIPTION

EC 214 and EC 215 are digital bench-top multirange conductivity meters designed for simplicity of use in measuring electrical conductivity in liquids.

Four ranges of conductivity measurements are provided to cover every application from deionized water to brine.

The calibration is made simple through the easy-to-operate front panel knob for slope adjustment.

The conductivity probe does not require re-calibration when switching from one range to another.


The probe is made of glass with 4 platinum rings that are corrosion resistant. It is also suitable for measuring conductivity of liquids in small sample sizes, and comes with a 1 m (3.3') cable.

The temperature effect can be compensated through a knob on the front panel of EC 214.

With EC 215, the 4-ring probe has a built-in temperature sensor that automatically compensates for temperature changes in the liquid tested. The temperature coefficient can be adjusted from 0 to 2.5% through a knob on the front panel.

It is also available a model with recorder output feature: EC 215R.

FUNCTIONAL DESCRIPTION EC 214




- 1) Liquid Crystal Display (LCD)
- 2) Conductivity calibration knob
- 3) 199.9 µS key, to select the range 0.0 to 199.9µS/cm
- 4) 19.99 mS key, to select the range 0.00 to 19.99 mS/cm
- 5) 199.9 mS key, to select the range 0.0 to 199.9mS/cm
- 6) 1999 µS key, to select the range 0 to 1999µS/cm
- 7) COND/□C key, to select conductivity reading or temperature setting for manual compensation
- 8) Temperature setting knob for manual compensation
- 9) ON/OFF switch
- 10) Power adapter socket
- 11) Probe connector

SPECIFICATIONS EC 214

Range	0.0 to 199.9 µS/cm 0 to 1999 µS/cm 0.00 to 19.99 mS/cm 0.0 to 199.9 mS/cm
Resolution	0.1 µS/cm 1 µS/cm 0.01 mS/cm 0.1 mS/cm
Accuracy (@20°C / 68°F)	±1% Full Scale (excluding probe error)
Typical EMC Deviation	±1% Full Scale
Calibration	Manual single point through front knob
Temperature Compensation	Manual from 0 to 50°C with B fixed at 2%/°C
Probe	HI 76300 platinum 4-ring conductivity probe with 1m (3.3') cable (included)
Power Supply	12 VDC adapter (HI710005 or HI710006 included)
Environment	0 to 50°C (32 to 122°F); 0-95% RH non-condensing
Dimensions	240x182x74 mm (9.4x7.1x2.9")
Weight	1.0 Kg (2.3 lb.)

FUNCTIONAL DESCRIPTION EC 215



- 1) Liquid Crystal Display (LCD)
- 2) Conductivity calibration knob
- 3) 199.9 μS key, to select the range 0.0 to 199.9 μS/cm
- 4) 19.99 mS key, to select the range 0.00 to 19.99 mS/cm
- 5) 199.9 mS key, to select the range 0.0 to 199.9 mS/cm
- 6) 1999 μS key, to select the range 0 to 1999 μS/cm
- 7) Temperature Coefficient setting knob
- 8) ON/OFF switch
- 9) Power adapter socket
- 10) Probe connector
- 11) Analog Output (EC 215R only)

SPECIFICATIONS EC 215

Range	0.0 to 199.9 μS/cm 0 to 1999 μS/cm 0.00 to 19.99 mS/cm 0.0 to 199.9 mS/cm
Resolution	0.1 μS/cm 1 μS/cm 0.01 mS/cm 0.1 mS/cm
Accuracy (@20°C / 68°F)	±1% Full Scale (excluding probe error)
Typical EMC Deviation	±1% Full Scale
Calibration	Manual at single point, through front knob
Temperature Compensation	Automatic from 0 to 50 °C with adjustable β from 0 to 2.5%/°C
Probe (included)	HI 76303 platinum 4-ring conductivity probe with temperature sensor and 1m (3.3') cable
Analog Output (EC 215R only)	0 to 5VCC not isolated analog output - 20mA max Range Output 0.0 to 199.9 μS/cm 0 to 5 VCC 0 to 1999 μS/cm 0 to 5 VCC 0.00 to 19.99 mS/cm 0 to 5 VCC 0.0 to 199.9 mS/cm 0 to 5 VCC
Analog Output Accuracy (EC 215R only)	0.1% of maximum displayed value
Analog Output Resolution (EC 215R only)	±2.5 mV
Power Supply	12 VDC power adapter (included)
Environment	0 to 50 °C (32 to 122 °F); 0-95% RH non-condensing
Dimensions	240x182x74 mm (9.4x7.1x2.9")
Weight	1.0 Kg (2.3 lb.)

OPERATIONAL GUIDE

Power connection

Plug the 12VDC adapter (HI710005 or HI710006) into the power supply socket (#10 on page 4; #9 on page 6).

Note: make sure the main line is protected by a fuse.

Probe connection


Connect the conductivity probe to the socket provided (#11 on page 4; #10 on page 6).

Note: Make sure that the instrument has been calibrated before taking conductivity measurements (see calibration paragraph).

Note: If possible, use plastic beakers to minimize EMC interferences.

TAKING CONDUCTIVITY MEASUREMENTS WITH EC 214

- switch the instrument on by pressing the ON/OFF key.
- Immerse the probe in the solution submerging the holes of the sleeve (0.5 cm below). Tap the probe lightly on the bottom of the beaker to remove any air bubbles trapped inside the sleeve. Take the temperature of the solution with a CHECKTEMP or a glass thermometer.
- Press the "COND/□C" key to select temperature setting.
- Adjust the "TEMPERATURE" knob to display the temperature of the solution on the LCD.



- Press the "COND/□C" key to select conductivity reading.
- Select the appropriate conductivity range.



8


Note: If the display shows "1", there is an over-range condition. Select the next higher range.

- Allow a few minutes for the reading to stabilize. The LCD will display the temperature compensated conductivity reading.



TAKING CONDUCTIVITY MEASUREMENTS WITH EC 215

- switch the instrument on by pressing the ON/OFF key.
- Immerse the probe in the solution submerging the holes of the sleeve (0.5 cm below). Tap the probe lightly on the bottom of the beaker to remove any air bubbles trapped inside the sleeve.
- Adjust the "TEMPERATURE COEFFICIENT" knob to the temperature coefficient value of the solution (see determination of the temperature coefficient of a solution paragraph).



- Select the appropriate conductivity range.



Note: If the display shows "1", there is an over-range condition. Select the next higher range.

- Allow a few minutes for the reading to stabilize. The LCD will display the temperature compensated conductivity reading.



9

CALIBRATION

The instrument should be re-calibrated at least once a month, or when the probe is changed.

INITIAL PREPARATION

If you are measuring in the mS ranges, calibrate the meter using HI7030 (or HI8030) conductivity solution (12.88 mS @ 25°C) or HI7034 (or HI8034) conductivity solution (80 mS @ 25°C). For the µS ranges, use HI7031 (or HI8031) conductivity solution (1413 µS @ 25°C) when calibrating in the range from 0 to 1999 µS or HI7033 (or HI8033) conductivity solution (84 µS @ 25°C) when calibrating in the range from 0 to 199.9 µS.


Choose a solution with a conductivity value close to the solution to be measured.

Rinse the probe thoroughly in distilled water. This is to minimize contamination of the calibration solution and secure higher accuracy.

When possible, use plastic beakers to minimize any EMC interferences.

PROCEDURE FOR EC 214:

- Pour a small quantity of the conductivity solution, e.g. HI7030 or HI8030 (12.88 mS/cm @ 25°C), into a plastic beaker.
- Immerse the probe in the solution submerging the holes of the sleeve (0.5 cm below). Tap the probe lightly on the bottom of the beaker to remove any air bubbles trapped inside the sleeve. Take the temperature of the solution with a CHECKTEMP or a glass thermometer.
- Press the "COND/°C" key to select temperature setting.
- Adjust the "TEMPERATURE" knob to display the temperature of the solution on the LCD.



10


- Press the "COND/°C" key to select conductivity reading and select the appropriate range:

"199.9 µS" for HI7033/HI8033,
 "1999 µS" for HI7031/HI8031,
 "19.99 mS" for HI7030/HI8030,
 "199.9 mS" for HI7034/HI8034.



- Note: If the display shows "1", there is an over-range condition. Select the next higher range.


- Allow a few minutes for the reading to stabilize and adjust the "CALIBRATION" knob to read (on the LCD) the value of the buffer solution at 25°C (77°F), e.g. 12.88 mS/cm.



- All subsequent measurements will be referenced to 25°C (77°F).
- To reference the measurements to 20°C (68°F), adjust the "CALIBRATION" knob to read (on the LCD) the value of the buffer solution at 20°C (68°F), e.g. 11.67 mS/cm. See conductivity vs. temperature chart.
- Calibration is now complete and the instrument is ready for use.

PROCEDURE FOR EC 215:

- Pour a small quantity of the conductivity solution, e.g. HI7030 or HI8030 (12.88 mS/cm @ 25°C), into a plastic beaker.
- Immerse the probe in the solution submerging the holes of the sleeve (0.5 cm below). Tap the probe lightly on the bottom of the beaker to remove any air bubbles trapped inside the sleeve.
- Adjust the "TEMPERATURE COEFFICIENT" knob to 2%/°C.




11

- Select the appropriate range:
 - "199.9 μ S" for HI7033/HI8033,
 - "1999 μ S" for HI7031/HI8031,
 - "19.99 mS" for HI7030/HI8030,
 - "199.9 mS" for HI7034/HI8034.

Note: If the display shows "1", there is an over-range condition. Select the next higher range.

- Allow a few minutes for the reading to stabilize and adjust the "CALIBRATION" knob to read (on the LCD) the value of the buffer solution at 25°C (77°F), e.g. 12.88 mS/cm.



- All subsequent measurements will be referenced to 25°C (77°F).
- To reference the measurements to 20°C (68°F), adjust the "CALIBRATION" knob to read (on the LCD) the value of the buffer solution at 20°C (68°F), e.g. 11.67 mS/cm. See conductivity vs. temperature chart.
- Calibration is now complete and the instrument is ready for use.

CONDUCTIVITY VERSUS TEMPERATURE CHART

The conductivity of an aqueous solution is the measure of its ability to carry an electrical current by means of ionic motion.

The conductivity invariably increases with increasing temperature.

It is affected by the type and number of ions in the solution and by the viscosity of the solution itself. Both parameters are temperature dependent. The dependency of conductivity on temperature is expressed as a relative change per degree Celsius at a particular temperature, commonly as percent per °C.

Since a small difference in temperature causes a large change in conductivity, the readings are usually normalized at 25°C.

°C	°F	HI 7030 HI 8030 (μ S/cm)	HI 7031 HI 8031 (μ S/cm)	HI 7033 HI 8033 (μ S/cm)	HI 7034 HI 8034 (μ S/cm)	HI 7035 HI 8035 (μ S/cm)	HI 7039 HI 8039 (μ S/cm)
0	32	7150	776	64	48300	65400	2760
5	41	8220	896	65	53500	74100	3180
10	50	9330	1020	67	59600	83200	3615
15	59	10480	1147	68	65400	92500	4063
16	60.8	10720	1173	70	67200	94400	4155
17	62.6	10950	1199	71	68500	96300	4245
18	64.4	11190	1225	73	69800	98200	4337
19	66.2	11430	1251	74	71300	100200	4429
20	68	11670	1278	76	72400	102100	4523
21	69.8	11910	1305	78	74000	104000	4617
22	71.6	12150	1332	79	75200	105900	4711
23	73.4	12390	1359	81	76500	107900	4805
24	75.2	12640	1386	82	78300	109800	4902
25	77	12880	1413	84	80000	111800	5000
26	78.8	13130	1440	86	81300	113800	5096
27	80.6	13370	1467	87	83000	115700	5190
28	82.4	13620	1494	89	84900	117700	5286
29	84.2	13870	1521	90	86300	119700	5383
30	86	14120	1548	92	88200	121800	5479
31	87.8	14370	1575	94	90000	123900	5575

EC 214 manually compensates for temperature differences with a fixed β at 2%.

EC 215 automatically compensates for temperature differences with a built-in circuitry. A knob is provided to adjust the temperature coefficient manually from 0% (without compensation) to 2.5% per degree Celsius.

DETERMINATION OF THE TEMPERATURE COEFFICIENT OF A SOLUTION (EC 215 & EC 215R)

- Immerse the probe into a sample of the solution and adjust the "TEMPERATURE COEFFICIENT" knob to 0% (i.e. no compensation).
- Condition the sample and probe at 25°C and note the conductivity reading C₂₅.
- Condition the sample and probe to a temperature t°C which is approximately 5°C to 10°C different from 25°C and note the conductivity reading C_t.
- The temperature coefficient β of the solution is calculated with the formula:

$$\beta = 100 \times \frac{(C_t - C_{25})}{(t - 25) \times C_{25}}$$

The above procedure is suitable for determining the temperature coefficient in the laboratory where the temperature of the solution can be determined and controlled. If this is not possible, e.g. during on-site measurements, the following procedure should be used:


- Immerse the probe into the test solution and turn the "TEMPERATURE COEFFICIENT" knob to 0% (no compensation).
- Allow the conductivity reading to stabilize (the reading should not change by more than ±0.2 mS within 1 minute) and record the value, C.
- Repeat the procedure with the temperature of the solution changed by more than 10°C. Wait for the conductivity reading to stabilize.
- Adjust the "TEMPERATURE COEFFICIENT" knob until the display reads the value C as recorded earlier.
- The value indicated by the knob is the temperature coefficient of the solution.

PROBE MAINTENANCE

Rinse the probe with tap water after every series of measurements. If a more thorough cleaning is required, remove the sleeve and clean the probe with a cloth or a nonabrasive detergent.

After cleaning the probe, re-calibrate the instrument.

The four ring platinum probe body is in glass. For this reason great care while handling the probe must be taken.



ACCESSORIES

CONDUCTIVITY BUFFER SOLUTIONS:

HI 7030L	12880 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460mL
HI 7030M	12880 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 230mL
HI 7031L	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460mL
HI 7031M	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 230mL
HI 7033L	84 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL
HI 7033M	84 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 230 mL
HI 7034L	80000 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460mL
HI 7034M	80000 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 230mL
HI 7035L	111800 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460mL
HI 7035M	111800 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 230mL
HI 7039L	5000 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460mL
HI 7039M	5000 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 230mL

CONDUCTIVITY BUFFER SOLUTIONS IN FDA APPROVED

BOTTLES:

HI 8030L	12880 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL
HI 8031L	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL
HI 8033L	84 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL
HI 8034L	80000 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL
HI 8035L	111800 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL
HI 8039L	5000 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$), 460 mL

ELECTRODE CLEANING SOLUTIONS

HI 7061M	General Cleaning Sol., 230 mL
HI 7061L	General Cleaning Sol., 460 mL

ELECTRODE CLEANING SOLUTIONS IN FDA APPROVED

BOTTLES


HI 8061M	General Cleaning Sol., 230 mL
HI 8061L	General Cleaning Sol., 460 mL

CONDUCTIVITY PROBES:

HI 76300	platinum 4-ring conductivity probe with 1m (3.3') cable (for EC214).
HI 76303	platinum 4-ring conductivity probe with temperature sensor and 1m (3.3') cable (for EC215).

OTHER ACCESSORIES:

CHECKTEMPC	Electronic thermometer (range: -50.0 to 150.0°C)
HI710005	Voltage adapter from 115VAC to 12VDC
HI710006	Voltage adapter from 230VAC to 12VDC



HI 7400011	Base for bench meters
HI 740036	100mL plastic beakers (6 pcs)
HI 740034	Cap for 100mL beakers (6 pcs)
HI 740134	Dust cover
HI 76405	Electrode holder

WARRANTY

All Hanna Instruments meters are guaranteed for two years against defects in workmanship and materials when used for their intended purpose and maintained according to instructions. The electrodes and the probes are guaranteed for a period of six months. This warranty is limited to repair or replacement free of charge.

Damage due to accident, misuse, tampering or lack of prescribed maintenance are not covered.

If service is required, contact the dealer from whom you purchased the instrument. If under warranty, report the model number, date of purchase, serial number and the nature of the failure. If the repair is not covered by the warranty, you will be notified of the charges incurred. If the instrument is to be returned to Hanna Instruments, first obtain a Returned Goods Authorization number from the Customer Service department and then send it with shipping costs prepaid. When shipping any instrument, make sure it is properly packaged for complete protection.

To validate your warranty, fill out and return the enclosed warranty card within 14 days from the date of purchase.

All rights are reserved. Reproduction in whole or in part is prohibited without the written consent of the copyright owner, Hanna Instruments Inc., 584 Park East Drive, Woonsocket, Rhode Island, 02895, USA.

Hanna Instruments reserves the right to modify the design, construction and appearance of its products without advance notice.

CE DECLARATION OF CONFORMITY



DECLARATION OF CONFORMITY

We

Hanna Instruments Italia Srl
via E.Fermi, 10
35030 Sarmeola di Rubano - PD
ITALY

herewith certify that the bench-top EC meters

EC 214 EC215

have been tested and found to be in compliance with EMC Directive 89/336/EEC and Low Voltage Directive 73/23/EEC according to the following applicable normatives:

EN 50082-1: Electromagnetic Compatibility - Generic Immunity Standard
IEC 801-2: Electrostatic Discharge
IEC 801-3: RF Radiated
IEC 801-4: Fast Transient

EN 50081-1: Electromagnetic Compatibility - Generic Emission Standard
EN 55022: Radiated, Class B

EN61010-1: Safety requirements for electrical equipment for measurement, control and laboratory use

Date of Issue: 02.02.2001


A. Marsilio - Technical Director
On behalf of
Hanna Instruments S.r.l.

Recommendations for Users

Before using these products, make sure that they are entirely suitable for the environment in which they are used. Operation of these instruments in residential area could cause unacceptable interferences to radio and TV equipments, requiring the operator to take all necessary steps to correct interferences.

The metal band at the end of the sensor is sensitive to electrostatic discharges. Avoid touching this metal band at all times.

During calibration of instruments, ESD wrist straps should be worn to avoid possible damage to the sensor by electrostatic discharge.

Any variation introduced by the user to the supplied equipment may degrade the instruments' EMC performance. To avoid electrical shock, do not use these instruments when voltages at the measurement surface exceed 24VAC or 60VDC.

Use plastic beakers to minimize any EMC interferences.

To avoid damages or burns, do not perform any measurement in microwave ovens.

MANE21SR2 05/01



<http://www.hannainst.com>