

Capsule Filling and Suppository & Pessary Moulds


Adelphi Hand Split Moulds

Hand made for over 50 years these moulds manufactured from brass and double nickel plated are accurate to within 10% based on Witepsol base. Standard moulds are conical in shape and range from 0.25 gm to 8 gm. Torpedo shaped moulds are available to order for sizes 1 gm to 4 gm only.

Part Numbers:

All conical moulds start 6005, all torpedo moulds start 6006, the last 3 numbers denote number of cavities and size and are written in table.

No. of Cavities	Size in grams 6005+ Part No. 6006+						No. of Cavities
		0.25	0.5	1	2	4	
6	000	001	002	003	004	005	6
12	100	101	102	103	104	105	12
24	200	201	202	203	204	205	24
36	300	301	302	303	304		36
50	400	401	402	403	404		50
70			502	503	504		70
100	600						100



The Aponorm Aponorm Fix 100 Part No: 6010000

This simple hand operated machine is designed to separate, fill and close hard gelatine capsules of size 0 only, 100 pieces at a time.

A similar machine exists for size 1 capsules.

Part No: 6010002

Specifications:

Materials: Nickel plated metal and plastic

Dimensions: Size 0: 225 x 235 x 80 mm

Size 1: 220 x 180 x 80 mm

Operational weight: Size 0: 1.6 kg, Size 1: 1.2 kg



Aponorm 60

This simple hand operated machine is designed to separate, fill and close hard gelatine capsules of any size 000-5, 60 pieces at a time by the use of interchangeable plates.

Specifications:

Materials: Nickel plated metal and plastic

Dimensions: 220 x 180 x 80 mm

Operational weight: 1.1 kg

Part No.	Item	Size of Capsule	Notes
6010001	Aponorm 60 base unit	Any	One size of plate included
6010100	Interchangeable plates	000	
6010101	Interchangeable plates	00	
6010102	Interchangeable plates	0	
6010103	Interchangeable plates	1	
6010104	Interchangeable plates	2	
6010105	Interchangeable plates	3	
6010106	Interchangeable plates	4	
6015100	Capsule orientator	0	For 30 capsules each time
6015101	Capsule orientator	1	For 30 capsules each time

Capsule Filling

The Adelphi Profill Filling System

Complies with US & International Pharmacopoeias

The Adelphi Profill Capsule Filling machine fills up to 100 capsules at a time of all types, and is supplied with one set of change parts to fill one size of capsule. Additional change parts are available to fill other capsule sizes. Fill weights within +/- 5-7% of average are possible, based on using pre-measured powder for each batch and using the tamping device supplied with each machine. Each machine is supplied with orientator, powder spreader, 50 pin tamper, powder tray and video CD instructions.

The machines are easy to clean and made to GMP standards, all powder contact parts being made from stainless steel or anodised aluminium. Double-blind capsules, tablets and caplets can be made speedily using the patent design and optional inserters.

Heat resistant parts allow autoclaving. All but 2 parts are dishwasher-safe. All parts can be dried with hot air.

All Adelphi Profill Capsule fillers are supplied with orientators and it is therefore not possible to interchange all sizes on one basic machine. However it is possible to fill all sizes except 000 if other change parts and conversion kits are bought.

The table below explains, by reference to part number, what is required to be able to fill other sizes of capsules.

All part numbers to be prefixed 6011



Adelphi Profill Capsule Filling Machines

Part No.	Size of capsule	Part No.	Description
6011000	- Size 000	6011007	- Size 4
6011001	- Size 00	6011008	- Size 5
6011002	- Size 0	6011009	- Size 00el
6011003	- Size 1L	6011010	- Size 0el
6011004	- Size 1S	6011011	- Size AA
6011005	- Size 2	6011012	- Size AAel
6011006	- Size 3		

Change parts for Adelphi Capsule filling machines

Part No.	Change of parts	Part No.	Conversion kit
6011200	for Size 000	6011300	Size 00-1
6011201	for Size 00	6011301	Size 1-3
6011202	for Size 0	6011302	Size 4-5
6011203	for Size 1		
6011205	for Size 2		
6011206	for Size 3		
6011207	for Size 4		
6011208	for Size 5		

Machine Size	000	00	0	1L	1S	2	3	4	5
To fill 000		200	200	200	200	200	200	200	200
To fill 00	001		201	201	201 & 300	201 & 300	201 & 300	201 & 300	201 & 300
To fill 0	002	202		202	202 & 300	202 & 300	202 & 300	202 & 300	202 & 300
To fill 1L/S	003	203	203			203	203	203 & 301	203 & 301
To fill 2	005	205 & 301	205 & 301	205 & 301	205		205	205 & 10	205 & 10
To fill 3	006	206 & 301	206 & 301	206 & 301	206	206		206 & 10	206 & 10
To fill 4	007	207 & 302	207 & 302	207 & 302	207 & 302	207 & 302	207 & 302		207
To fill 5	008	208 & 302	208 & 302	208 & 302	208 & 302	208 & 302	208 & 302	208	

Capsule sizes that can be filled with Adelphi Capsule Filling System

Machine Size	000	00	0	1L	1S	2	3	4	5
With additional change parts	000	00,0,1	00,0,1	00,0,1	1,2,3	1,2,3	1,2,3	4,5	4,5
With additional change parts and conversion kit	N/A	Any size capsule from 000-5							