

Specifications - Volumetric KF Titrator Base EVA V1

Technique	Upol DC, amperometric
User rights management	Fully configurable
Measurement Type	Karl Fischer (Volumetric)
Data Transfer	Automatic sample weight transfer SmartChemicals Transfer data with barcode reader
Burette Dispensing Resolution	120,000 steps
Burette Dispensing Accuracy and Precision	According to DIN/ISO 8655-3
Sample Type	Liquids Soluble solids
Limit of Quantification (LOQ)	0.2 mg water per sample (RSD = 1%)
Sample Handling Technique	Single drying oven possibility
Product Line	Titration EVA