Application Range Ion-selective Measurements ● Recommended by WTW ○ Suitable		
Application Range	pH/ION 7320	pH/ION 340i/3400i*, Multi 350i/3500i*
Occasional, simple ISE measurement	0	•
Routine and standard measurement	•	0
Advanced methods and procedures	•	_
see page	49	22, 51

Laboratory ISE Benchtop Meters

Reliable documentation of ISE measurements...

... with the inoLab® pH/ION 7320

The new inoLab® pH/ION 7320 is ideal for precision measurements and automatic documentation acc. to GLP/ AQA in quality laboratories for all branches. Optional available with built-in printer.

inoLab® pH/ION 7320**

- 2-channel instrument for simultaneous measuring of pH, ISE or Redox
- Data transfer via USB interface
- For fast data transfer in .csv format or via the optional integrated printer
- CMC function for monitoring the measuring range for pH and ISE measurements
- * North American version
- ** available in Q4/2012

